

## HemaPro™

### Ex Vivo Reagent

- Ex Vivo Hematopoietic Stem Cell Expansion Solution
- Clonogenic Medium for the Culturing of Human Progenitor Cells
- Serum Free Solution Contains
  - No erythropoietin
  - No recombinant growth factors
  - No human serum or fetal bovine serum
  - No mitogenic contaminants
- Effective in Studying Stimulatory Factors

### Specifications

- Classification – Defined
- Stock Concentration – 50X
- Supplement to liquid basal media
- Recommended Storage – 2-8 °C, dark
- Stability – Minimum 6 months at 50X
- Stability – 45 days when added to basal medium

**For Laboratory use. Not For injection**